

For Research Use Only

# ARID1A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83733-2-PBS

## Basic Information

Catalog Number:

83733-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30749

GenBank Accession Number:

NM\_006015

GeneID (NCBI):

8289

UNIPROT ID:

O14497

Full Name:

AT rich interactive domain 1A (SWI-like)

Calculated MW:

242 kDa

Purification Method:

Protein A purification

CloneNo.:

240690G2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00714-2, ARID1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83733-3-PBS. Detection antibody: 83733-2-PBS. Standard: Ag30749. Range: 1.25 -80 ng/mL.