

For Research Use Only

USP4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83734-3-PBS

Basic Information

Catalog Number: 83734-3-PBS	GenBank Accession Number: BC125131	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7375	CloneNo.: 240797B6
Source: Rabbit	UNIPROT ID: Q13107	
Isotype: IgG	Full Name: ubiquitin specific peptidase 4 (proto-oncogene)	
Immunogen Catalog Number: AG28297	Calculated MW: 963 aa, 109 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00717-2, USP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83734-3-PBS. Detection antibody: 83734-1-PBS. Standard: Ag28297. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00717-3, USP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83734-3-PBS. Detection antibody: 83734-4-PBS. Standard: Ag28297. Range: 0.625-80 ng/mL