

For Research Use Only

USP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83734-4-PBS



Basic Information

Catalog Number:

83734-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28297

GenBank Accession Number:

BC125131

GeneID (NCBI):

7375

UNIPROT ID:

Q13107

Full Name:

ubiquitin specific peptidase 4 (proto-oncogene)

Calculated MW:

963 aa, 109 kDa

Purification Method:

Protein A purification

CloneNo.:

240797G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00717-1, USP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83734-2-PBS. Detection antibody: 83734-4-PBS. Standard: Ag28297. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00717-3, USP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83734-3-PBS. Detection antibody: 83734-4-PBS. Standard: Ag28297. Range: 0.625-80 ng/mL



Bi-layer interferometry (BLI) kinetic assay of 83734-4-PBS against Human USP4 was performed. The affinity constant is 16.1 nM.