

For Research Use Only

CBR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83736-2-PBS



Basic Information

Catalog Number:

83736-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34830

GenBank Accession Number:

BC002511

GeneID (NCBI):

873

UNIPROT ID:

P16152

Full Name:

carbonyl reductase 1

Calculated MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

240707F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00724-1, CBR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83736-2-PBS. Detection antibody: 83736-1-PBS. Standard: Ag34830. Range: 1.56-100 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83736-2-PBS against Human CBR1 was performed. The affinity constant is 76.2 pM.