

For Research Use Only

CLEC7A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83739-2-PBS

Basic Information

Catalog Number: 83739-2-PBS	GenBank Accession Number: BC013385	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64581	CloneNo.: 240765G3
Source: Rabbit	UNIPROT ID: Q9BXN2	
Isotype: IgG	Full Name: C-type lectin domain family 7, member A	
Immunogen Catalog Number: AG18924	Calculated MW: 247 aa, 28 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00689-1, Dectin-1/CLEC7A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83739-2-PBS. Detection antibody: 83739-1-PBS. Standard: Ag18924. Range: 0.78-100 ng/mL.