

For Research Use Only

RAB37 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83740-2-PBS

Basic Information

Catalog Number: 83740-2-PBS	GenBank Accession Number: BC016615	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 326624	CloneNo.: 240838D9
Source: Rabbit	UNIPROT ID: Q96AX2	
Isotype: IgG	Full Name: RAB37, member RAS oncogene family	
Immunogen Catalog Number: AG3717	Calculated MW: 223 aa, 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00690-1, RAB37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83740-3-PBS. Detection antibody: 83740-2-PBS. Standard: Ag3717. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00690-2, RAB37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83740-1-PBS. Detection antibody: 83740-2-PBS. Standard: Ag3717. Range: 0.625-80 ng/mL