

For Research Use Only

Mouse RAGE Recombinant antibody, PBS Only (Detector)



Catalog Number: 83742-1-PBS

Basic Information

Catalog Number:

83742-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007425.3

GeneID (NCBI):

11596

UNIPROT ID:

Q62151-1

Full Name:

advanced glycosylation end product-specific receptor

Calculated MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

240774A12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00693-1, MOUSE RAGE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83742-2-PBS. Detection antibody: 83742-1-PBS. Standard: Eg0873. Range: 0.78-100 ng/mL.



Cytometric bead array standard curve of MP00693-2, MOUSE RAGE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83742-4-PBS. Detection antibody: 83742-1-PBS. Standard: Eg0873. Range: 0.78-100 ng/mL.