

For Research Use Only

# MSLN Recombinant antibody, PBS Only (Detector)



Catalog Number: 83744-2-PBS

## Basic Information

Catalog Number:

83744-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005823

GeneID (NCBI):

10232

UNIPROT ID:

Q13421

Full Name:

mesothelin

Calculated MW:

69KD

Purification Method:

Protein A purification

CloneNo.:

240795G5

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00697-2, MSLN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83744-5-PBS. Detection antibody: 83744-2-PBS. Standard: Eg31658. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00697-3, MSLN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83744-3-PBS. Detection antibody: 83744-2-PBS. Standard: Eg31658. Range: 0.625-40 ng/mL