

For Research Use Only

Rat Cxcl1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83745-1-PBS

Basic Information

Catalog Number: 83745-1-PBS	GenBank Accession Number: NM_030845.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81503	CloneNo.: 240799B2
Source: Rabbit	UNIPROT ID: P14095	
Isotype: IgG	Full Name: chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)	
	Calculated MW: 10 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00698-1, RAT CXCL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83745-1-PBS. Detection antibody: 83745-2-PBS. Standard: Eg0948. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00698-2, RAT CXCL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83745-1-PBS. Detection antibody: 83745-3-PBS. Standard: Eg0948. Range: 0.313-40 ng/mL.