

For Research Use Only

Rat Cxcl1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83745-2-PBS

Basic Information

Catalog Number:

83745-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_030845.1

GeneID (NCBI):

81503

UNIPROT ID:

P14095

Full Name:

chemokine (C-X-C motif) ligand 1
(melanoma growth stimulating
activity, alpha)

Calculated MW:

10 kDa

Purification Method:

Protein A purification

CloneNo.:

240799D7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00698-1, RAT CXCL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83745-1-PBS. Detection antibody: 83745-2-PBS. Standard: Eg0948. Range: 0.313-40 ng/mL.