

For Research Use Only

SLC35A1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83746-1-PBS

Basic Information

Catalog Number:

83746-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9524

GenBank Accession Number:

BC017807

GeneID (NCBI):

10559

UNIPROT ID:

P78382

Full Name:

solute carrier family 35 (CMP-sialic acid transporter), member A1

Calculated MW:

337 aa, 37 kDa

Purification Method:

Protein A purification

CloneNo.:

240824F6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00696-1, SLC35A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83746-2-PBS. Detection antibody: 83746-1-PBS. Standard: Ag9524. Range: 0.78-100 ng/mL.