

For Research Use Only

IL6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83747-2-PBS

Basic Information

Catalog Number: 83747-2-PBS	GenBank Accession Number: BC015511	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3569	CloneNo.: 240670A4
Source: Rabbit	ENSEMBL Gene ID: ENSG00000136244	
Isotype: IgG	UNIPROT ID: P05231	
	Full Name: interleukin 6 (interferon, beta 2)	
	Calculated MW: 212 aa, 24 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00687-1, IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83747-2-PBS. Detection antibody: 83747-4-PBS. Standard: Eg0117. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00687-2, IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83747-2-PBS. Detection antibody: 83747-1-PBS. Standard: Eg0117. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00687-3, IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83747-2-PBS. Detection antibody: 83747-3-PBS. Standard: Eg0117. Range: 0.625-80 ng/mL