

For Research Use Only

LIN28 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83748-3-PBS

Basic Information

Catalog Number: 83748-3-PBS	GenBank Accession Number: BC028566	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79727	CloneNo.: 240683C4
Source: Rabbit	UNIPROT ID: Q9H9Z2	
Isotype: IgG	Full Name: lin-28 homolog (C. elegans)	
Immunogen Catalog Number: AG2312	Calculated MW: 209 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00699-2, LIN28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83748-3-PBS. Detection antibody: 83748-1-PBS. Standard: Ag2312. Range: 0.78-100 ng/mL