

For Research Use Only

# CKMT1A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83750-3-PBS

## Basic Information

|                                            |                                                        |                                                       |
|--------------------------------------------|--------------------------------------------------------|-------------------------------------------------------|
| <b>Catalog Number:</b><br>83750-3-PBS      | <b>GenBank Accession Number:</b><br>BC001926           | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                    | <b>GeneID (NCBI):</b><br>548596                        | <b>CloneNo.:</b><br>240825E11                         |
| <b>Source:</b><br>Rabbit                   | <b>UNIPROT ID:</b><br>P12532                           |                                                       |
| <b>Isotype:</b><br>IgG                     | <b>Full Name:</b><br>creatine kinase, mitochondrial 1A |                                                       |
| <b>Immunogen Catalog Number:</b><br>AG7583 | <b>Calculated MW:</b><br>47 kDa                        |                                                       |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00706-2, CKMT1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83750-1-PBS. Detection antibody: 83750-3-PBS. Standard: Ag7583. Range: 0.313-40 ng/mL.