

For Research Use Only

# GBP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83751-2-PBS

## Basic Information

Catalog Number: 83751-2-PBS	GenBank Accession Number: BC002666	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2633	CloneNo.: 240634G7
Source: Rabbit	UNIPROT ID: P32455	
Isotype: IgG	Full Name: GTP binding protein 1	
Immunogen Catalog Number: AG7562	Calculated MW: 68 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00707-1, GBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83751-2-PBS. Detection antibody: 83751-1-PBS. Standard: Ag7562. Range: 0.78-100 ng/mL