

For Research Use Only

GBP1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83751-5-PBS

Basic Information

Catalog Number: 83751-5-PBS	GenBank Accession Number: BC002666	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2633	CloneNo.: 240634F7
Source: Rabbit	UNIPROT ID: P32455	
Isotype: IgG	Full Name: GTP binding protein 1	
Immunogen Catalog Number: AG7562	Calculated MW: 68 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00707-2, GBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83751-3-PBS. Detection antibody: 83751-5-PBS. Standard: Ag7562. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00707-3, GBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83751-4-PBS. Detection antibody: 83751-5-PBS. Standard: Ag7562. Range: 0.78-100 ng/mL