

For Research Use Only

# GBP1 Recombinant antibody, PBS Only



Catalog Number: 83751-6-PBS

## Basic Information

Catalog Number:

83751-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7562

GenBank Accession Number:

BC002666

GeneID (NCBI):

2633

UNIPROT ID:

P32455

Full Name:

GTP binding protein 1

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

240634G8

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83751-6-RR (GBP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83751-6-PBS in a different storage buffer formulation.