

For Research Use Only

# COL1A2 Recombinant antibody, PBS Only (Detector)



Catalog Number:83752-1-PBS

## Basic Information

|  |  |   |
|--|--|---|
| <b>Catalog Number:</b><br>83752-1-PBS      | <b>GenBank Accession Number:</b><br>BC054498   | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                    | <b>GeneID (NCBI):</b><br>1278                  | <b>CloneNo.:</b><br>240665B7                          |
| <b>Source:</b><br>Rabbit                   | <b>UNIPROT ID:</b><br>P08123                   |   |
| <b>Isotype:</b><br>IgG                     | <b>Full Name:</b><br>collagen, type I, alpha 2 |   |
| <b>Immunogen Catalog Number:</b><br>AG6281 | <b>Calculated MW:</b><br>1366 aa, 130 kDa      |   |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00716-1, COL1A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83752-3-PBS. Detection antibody: 83752-1-PBS. Standard: Ag6281. Range: 1.56-100 ng/mL.