

For Research Use Only

# DCN Recombinant antibody, PBS Only (Capture)



Catalog Number: 83753-1-PBS

## Basic Information

Catalog Number: 83753-1-PBS	GenBank Accession Number: BC005322	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1634	CloneNo.: 240781E2
Source: Rabbit	UNIPROT ID: P07585	
Isotype: IgG	Full Name: decorin	
Immunogen Catalog Number: AG6275	Calculated MW: 359 aa, 40 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00719-2, Decorin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83753-1-PBS. Detection antibody: 83753-2-PBS. Standard: Ag6275. Range: 1.56-100 ng/mL