

For Research Use Only

DCN Recombinant antibody, PBS Only (Detector)



Catalog Number: 83753-3-PBS

Basic Information

Catalog Number: 83753-3-PBS	GenBank Accession Number: BC005322	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1634	CloneNo.: 240781G5
Source: Rabbit	UNIPROT ID: P07585	
Isotype: IgG	Full Name: decorin	
Immunogen Catalog Number: AG6275	Calculated MW: 359 aa, 40 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00719-1, Decorin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83753-2-PBS. Detection antibody: 83753-3-PBS. Standard: Ag6275. Range: 1.56-100 ng/mL