

For Research Use Only

# DCN Recombinant antibody, PBS Only



Catalog Number: 83753-4-PBS

## Basic Information

Catalog Number: 83753-4-PBS	GenBank Accession Number: BC005322	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1634	CloneNo.: 240781E3
Source: Rabbit	UNIPROT ID: P07585	
Isotype: IgG	Full Name: decorin	
Immunogen Catalog Number: AG6275	Calculated MW: 359 aa, 40 kDa	

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83753-4-RR (Decorin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83753-4-PBS in a different storage buffer formulation.