

For Research Use Only

BCL9 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83760-2-PBS

Basic Information

Catalog Number: 83760-2-PBS	GenBank Accession Number: BC116451	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 607	CloneNo.: 240682G8
Source: Rabbit	UNIPROT ID: O00512	
Isotype: IgG	Full Name: B-cell CLL/lymphoma 9	
Immunogen Catalog Number: AG19158	Calculated MW: 1426 aa, 149 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00739-1, BCL9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83760-1-PBS. Detection antibody: 83760-2-PBS. Standard: Ag19158. Range: 0.313-40 ng/mL