

For Research Use Only

BCL9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83760-6-PBS



Basic Information

Catalog Number: 83760-6-PBS	GenBank Accession Number: BC116451	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 607	CloneNo.: 240682E12
Source: Rabbit	UNIPROT ID: O00512	
Isotype: IgG	Full Name: B-cell CLL/lymphoma 9	
Immunogen Catalog Number: AG19158	Calculated MW: 1426 aa, 149 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00739-4, Human BCL9 Recombinant Matched Antibody Pair - PBS only. 83760-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19158. 83760-5-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean BCL9 concentration was determined to be 18.2 ng/mL in COLO 320 cell extract based on a 2.0 mg/mL extract load.