

For Research Use Only

# PAK4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83761-2-PBS

## Basic Information

<b>Catalog Number:</b> 83761-2-PBS	<b>GenBank Accession Number:</b> BC025282	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 10298	<b>CloneNo.:</b> 240686E1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> O96013	
<b>Isotype:</b> IgG	<b>Full Name:</b> p21 protein (Cdc42/Rac)-activated kinase 4	
<b>Immunogen Catalog Number:</b> AG5981	<b>Calculated MW:</b> 64 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00743-1, PAK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83761-3-PBS. Detection antibody: 83761-2-PBS. Standard: Ag5981. Range: 1.25-80 ng/mL