

For Research Use Only

PAK4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83761-7-PBS



Basic Information

Catalog Number: 83761-7-PBS	GenBank Accession Number: BC025282	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10298	CloneNo.: 240686E12
Source: Rabbit	UNIPROT ID: O96013	
Isotype: IgG	Full Name: p21 protein (Cdc42/Rac)-activated kinase 4	
Immunogen Catalog Number: AG5981	Calculated MW: 64 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00743-3, Human PAK4 Recombinant Matched Antibody Pair - PBS only. 83761-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5981. 83761-6-PBS was HRP conjugated as the detection antibody. Range: 0.469-30 ng/mL

HeLa cell lysate was measured. The human PAK4 concentration of detected samples was determined to be 50.68 ng/mL (based on a 1.20 mg/mL extract load).