

For Research Use Only

# Granzyme B Recombinant antibody, PBS Only (Capture)

Catalog Number:83763-2-PBS



## Basic Information

Catalog Number:

83763-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC030195

GeneID (NCBI):

3002

UNIPROT ID:

P10144

Full Name:

granzyme B (granzyme 2, cytotoxic T-lymphocyte-associated serine esterase 1)

Calculated MW:

247 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

240761F3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00761-2, Granzyme B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83763-2-PBS. Detection antibody: 83763-4-PBS. Standard: Eg0592. Range: 0.625-80 ng/mL.