

For Research Use Only

HIF-1 alpha/HIF1A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83764-1-PBS

Basic Information

Catalog Number: 83764-1-PBS	GenBank Accession Number: NM_001530.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3091	CloneNo.: 240767E7
Source: Rabbit	UNIPROT ID: Q16665-1	
Isotype: IgG	Full Name: hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor)	
	Calculated MW: 93 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00762-1, HIF-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83764-2-PBS. Detection antibody: 83764-1-PBS. Standard: Eg0939. Range: 0.625-80 ng/mL.