

For Research Use Only

IL33 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83765-1-PBS

Basic Information

| | | |
|--------------------------------|---------------------------------------|--|
| Catalog Number: 83765-1-PBS | GenBank Accession Number: BC047085 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 90865 | CloneNo.: 240719C7 |
| Source: Rabbit | UNIPROT ID: O95760 | |
| Isotype: IgG | Full Name: interleukin 33 | |
| | Calculated MW: 30 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00763-2, IL-33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83765-2-PBS. Detection antibody: 83765-1-PBS. Standard: Eg0313. Range: 6.25-100 ng/mL