

For Research Use Only

# IL33 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83765-3-PBS

## Basic Information

|                                |                                       |  |
|--------------------------------|---------------------------------------|--|
| Catalog Number:<br>83765-3-PBS | GenBank Accession Number:<br>BC047085 | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml               | GeneID (NCBI):<br>90865               | CloneNo.:<br>240719F5                          |
| Source:<br>Rabbit              | UNIPROT ID:<br>O95760                 |  |
| Isotype:<br>IgG                | Full Name:<br>interleukin 33          |  |
|                                | Calculated MW:<br>30 kDa              |  |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00763-1, IL-33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83765-2-PBS. Detection antibody: 83765-3-PBS. Standard: Eg0313. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP00763-3, IL-33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83765-4-PBS. Detection antibody: 83765-3-PBS. Standard: Eg0313. Range: 6.25-100 ng/mL.