

For Research Use Only

# VE-cadherin Recombinant antibody, PBS Only (Capture)



Catalog Number:83766-1-PBS

## Basic Information

Catalog Number:

83766-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001795.5

GeneID (NCBI):

1003

UNIPROT ID:

P33151-1

Full Name:

cadherin 5, type 2 (vascular endothelium)

Calculated MW:

88kDa

Purification Method:

Protein A purification

CloneNo.:

240755A6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00752-3, VE-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83766-1-PBS. Detection antibody: 83766-4-PBS. Standard: Eg1200. Range: 0.78-100 ng/mL