

For Research Use Only

VE-cadherin Recombinant antibody, PBS Only (Detector)



Catalog Number:83766-3-PBS

Basic Information

Catalog Number:

83766-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001795.5

GeneID (NCBI):

1003

UNIPROT ID:

P33151-1

Full Name:

cadherin 5, type 2 (vascular endothelium)

Calculated MW:

88kDa

Purification Method:

Protein A purification

CloneNo.:

240755C7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00752-2, VE-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83766-2-PBS. Detection antibody: 83766-3-PBS. Standard: Eg1200. Range: 0.78-100 ng/mL