

For Research Use Only

# TNC/Tenascin-C Recombinant antibody, PBS Only (Capture)

Catalog Number:83767-2-PBS



## Basic Information

Catalog Number: 83767-2-PBS	GenBank Accession Number: BC151843	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3371	CloneNo.: 240844D8
Source: Rabbit	UNIPROT ID: P24821	
Isotype: IgG	Full Name: tenascin C	
Immunogen Catalog Number: AG27122	Calculated MW: 2201 aa, 241 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00755-1, TNC/Tenascin-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83767-2-PBS. Detection antibody: 83767-4-PBS. Standard: Ag27122. Range: 0.78-100 ng/mL.