

For Research Use Only

Coagulation Factor III/Tissue Factor Recombinant antibody, PBS Only (Detector)

Catalog Number: 83776-3-PBS



Basic Information

Catalog Number: 83776-3-PBS	GenBank Accession Number: NM_001993.5	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2152	CloneNo.: 240784D6
Source: Rabbit	UNIPROT ID: P13726-1	
Isotype: IgG	Full Name: coagulation factor III (thromboplastin, tissue factor)	
	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00732-2, Coagulation Factor III/Tissue Factor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83776-1-PBS. Detection antibody: 83776-3-PBS. Standard: Eg1052. Range: 0.625-80 ng/mL.