

For Research Use Only

LILRA4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83777-1-PBS

Basic Information

Catalog Number: 83777-1-PBS	GenBank Accession Number: BC109198	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23547	CloneNo.: 240820B1
Source: Rabbit	UNIPROT ID: P59901	
Isotype: IgG	Full Name: leukocyte immunoglobulin-like receptor, subfamily A (with TM domain), member 4	
Immunogen Catalog Number: AG15721	Calculated MW: 499 aa, 55 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00735-1, LILRA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83777-2-PBS. Detection antibody: 83777-1-PBS. Standard: Ag15721. Range: 6.25-200 ng/mL