

For Research Use Only

CTLA4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83778-5-PBS

Basic Information

Catalog Number: 83778-5-PBS	GenBank Accession Number: BC070162	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1493	CloneNo.: 240793F7
Source: Rabbit	UNIPROT ID: P16410	
Isotype: IgG	Full Name: cytotoxic T-lymphocyte-associated protein 4	
	Calculated MW: 223 aa, 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00736-1, CTLA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83778-1-PBS. Detection antibody: 83778-5-PBS. Standard: Eg0972. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00736-2, CTLA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83778-4-PBS. Detection antibody: 83778-5-PBS. Standard: Eg0972. Range: 0.313-40 ng/mL