

For Research Use Only

ISLR Recombinant antibody, PBS Only (Detector)

Catalog Number:83782-4-PBS



Basic Information

Catalog Number: 83782-4-PBS	GenBank Accession Number: BC022478	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3671	CloneNo.: 240727C1
Source: Rabbit	UNIPROT ID: O14498	
Isotype: IgG	Full Name: immunoglobulin superfamily containing leucine-rich repeat	
Immunogen Catalog Number: AG19446	Calculated MW: 428 aa, 46 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00747-1, ISLR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83782-2-PBS. Detection antibody: 83782-4-PBS. Standard: Ag19446. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00747-2, ISLR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83782-3-PBS. Detection antibody: 83782-4-PBS. Standard: Ag19446. Range: 0.313-20 ng/mL



Bi-layer interferometry (BLI) kinetic assay of 83782-4-PBS against Human ISLR was performed. The affinity constant is 5.53 nM.