

For Research Use Only

# Islet 1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83784-3-PBS

## Basic Information

Catalog Number: 83784-3-PBS	GenBank Accession Number: BC031213	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3670	CloneNo.: 240833A9
Source: Rabbit	UNIPROT ID: P61371	
Isotype: IgG	Full Name: ISL LIM homeobox 1	
Immunogen Catalog Number: AG33541	Calculated MW: 39 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00771-2, Islet 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83784-1-PBS. Detection antibody: 83784-3-PBS. Standard: Ag33541. Range: 0.313- 40 ng/mL