

For Research Use Only

ATOX1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83785-3-PBS

Basic Information

Catalog Number: 83785-3-PBS	GenBank Accession Number: BC112248	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 475	CloneNo.: 240679E12
Source: Rabbit	UNIPROT ID: O00244	
Isotype: IgG	Full Name: ATX1 antioxidant protein 1 homolog (yeast)	
Immunogen Catalog Number: AG18460	Calculated MW: 68 aa, 8 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00767-2, ATOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83785-3-PBS. Detection antibody: 83785-1-PBS. Standard: Ag18460. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00767-3, ATOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83785-3-PBS. Detection antibody: 83785-5-PBS. Standard: Ag18460. Range: 0.625-80 ng/mL