

For Research Use Only

# HIF2 $\alpha$ /EPAS1 Recombinant antibody, PBS Only (Detector)

Catalog Number:83790-3-PBS



## Basic Information

Catalog Number: 83790-3-PBS	GenBank Accession Number: BC051338	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2034	CloneNo.: 240736F12
Source: Rabbit	UNIPROT ID: Q99814	
Isotype: IgG	Full Name: endothelial PAS domain protein 1	
Immunogen Catalog Number: AG15199	Calculated MW: 96 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00775-3, HIF2  $\alpha$  / EPAS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83790-4-PBS. Detection antibody: 83790-3-PBS. Standard: Ag15199. Range: 1.56-100 ng/mL.