

For Research Use Only

CEBPB Recombinant antibody, PBS Only (Detector)



Catalog Number:83791-1-PBS

Basic Information

Catalog Number: 83791-1-PBS	GenBank Accession Number: BC007538	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1051	CloneNo.: 240836G12
Source: Rabbit	UNIPROT ID: P17676	
Isotype: IgG	Full Name: CCAAT/enhancer binding protein (C/EBP), beta	
Immunogen Catalog Number: AG20073	Calculated MW: 345 aa, 36 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00737-1, CEBPB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83791-5-PBS. Detection antibody: 83791-1-PBS. Standard: Ag20073. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00737-2, CEBPB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83791-2-PBS. Detection antibody: 83791-1-PBS. Standard: Ag20073. Range: 0.313-40 ng/mL