

For Research Use Only

GPD1L Recombinant antibody, PBS Only (Capture)



Catalog Number:83792-1-PBS

Basic Information

Catalog Number: 83792-1-PBS	GenBank Accession Number: BC028726	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23171	CloneNo.: 240846B2
Source: Rabbit	UNIPROT ID: Q8N335	
Isotype: IgG	Full Name: glycerol-3-phosphate dehydrogenase 1-like	
Immunogen Catalog Number: AG10977	Calculated MW: 351 aa, 38 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00741-1, GPD1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83792-1-PBS. Detection antibody: 83792-2-PBS. Standard: Ag10977. Range: 1.56-100 ng/mL