

For Research Use Only

CDH17 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83793-4-PBS

Basic Information

Catalog Number: 83793-4-PBS	GenBank Accession Number: BC113464	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1015	CloneNo.: 240681F7
Source: Rabbit	UNIPROT ID: Q12864	
Isotype: IgG	Full Name: cadherin 17, LI cadherin (liver-intestine)	
Immunogen Catalog Number: AG19488	Calculated MW: 832 aa, 92 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00742-3, Cadherin-17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83793-4-PBS. Detection antibody: 83793-1-PBS. Standard: Ag19488. Range: 1.25-80 ng/mL.