

For Research Use Only

5 Lipoxygenase Recombinant antibody, PBS Only (Capture)

Catalog Number:83794-3-PBS



Basic Information

Catalog Number: 83794-3-PBS	GenBank Accession Number: BC130332	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 240	CloneNo.: 240740C1
Source: Rabbit	UNIPROT ID: P09917	
Isotype: IgG	Full Name: arachidonate 5-lipoxygenase	
Immunogen Catalog Number: AG19503	Calculated MW: 674 aa, 78 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00753-1, ALOX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83794-3-PBS. Detection antibody: 83794-2-PBS. Standard: Ag19503. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00753-2, ALOX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83794-3-PBS. Detection antibody: 83794-1-PBS. Standard: Ag19503. Range: 1.56-100 ng/mL