

For Research Use Only

# Mouse Il6 Recombinant antibody, PBS Only (Detector)



Catalog Number:83795-2-PBS

## Basic Information

Catalog Number:

83795-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_031168

GeneID (NCBI):

16193

UNIPROT ID:

P08505

Full Name:

interleukin 6

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

240701E6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00756-2, MOUSE IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83795-4-PBS. Detection antibody: 83795-2-PBS. Standard: Eg0973. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00756-3, MOUSE IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83795-5-PBS. Detection antibody: 83795-2-PBS. Standard: Eg0973. Range: 0.625-80 ng/mL.