

For Research Use Only

Mouse Cdkn2a/p16INK4a Recombinant antibody, PBS Only (Detector)

Catalog Number:83797-1-PBS



Basic Information

Catalog Number: 83797-1-PBS	GenBank Accession Number: NM_001040654	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12578	CloneNo.: 240935E2
Source: Rabbit	UNIPROT ID: Q64364	
Isotype: IgG	Full Name: cyclin-dependent kinase inhibitor 2A	
Immunogen Catalog Number: AG28927	Calculated MW: 18 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00764-2, MOUSE Cdkn2a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83797-2-PBS. Detection antibody: 83797-1-PBS. Standard: Ag28927. Range: 3.13-100 ng/mL.