

For Research Use Only

# Mouse Cdkn2a/p16INK4a Recombinant antibody, PBS Only (Capture)

Catalog Number:83797-4-PBS



## Basic Information

<b>Catalog Number:</b> 83797-4-PBS	<b>GenBank Accession Number:</b> NM_001040654	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 12578	<b>CloneNo.:</b> 240935E8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q64364	
<b>Isotype:</b> IgG	<b>Full Name:</b> cyclin-dependent kinase inhibitor 2A	
<b>Immunogen Catalog Number:</b> AG28927	<b>Calculated MW:</b> 18 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00764-3, MOUSE Cdkn2a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83797-4-PBS. Detection antibody: 83797-6-PBS. Standard: Ag28927. Range: 1.563-50 ng/mL.