

For Research Use Only

SYN2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83799-3-PBS

Basic Information

Catalog Number: 83799-3-PBS	GenBank Accession Number: BC051307	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6854	CloneNo.: 240852B1
Source: Rabbit	UNIPROT ID: Q92777	
Isotype: IgG	Full Name: synapsin II	
Immunogen Catalog Number: AG9855	Calculated MW: 582 aa, 63 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00727-1, SYN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83799-2-PBS. Detection antibody: 83799-3-PBS. Standard: Ag9855. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00727-3, SYN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83799-3-PBS. Detection antibody: 83799-1-PBS. Standard: Ag9855. Range: 0.781-50 ng/mL.