

For Research Use Only

IRX6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83800-2-PBS

Basic Information

Catalog Number: 83800-2-PBS	GenBank Accession Number: BC136573	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79190	CloneNo.: 240876F6
Source: Rabbit	UNIPROT ID: P78412	
Isotype: IgG	Full Name: iroquois homeobox 6	
Immunogen Catalog Number: AG19586	Calculated MW: 446 aa, 48 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00729-2, IRX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83800-2-PBS. Detection antibody: 83800-1-PBS. Standard: Ag19586. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00729-3, IRX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83800-2-PBS. Detection antibody: 83800-3-PBS. Standard: Ag19586. Range: 0.313-40 ng/mL