

For Research Use Only

IRX6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83800-3-PBS

Basic Information

| | | |
|---|--|---|
| Catalog Number: 83800-3-PBS | GenBank Accession Number: BC136573 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 79190 | CloneNo.: 240876F4 |
| Source: Rabbit | UNIPROT ID: P78412 | |
| Isotype: IgG | Full Name: iroquois homeobox 6 | |
| Immunogen Catalog Number: AG19586 | Calculated MW: 446 aa, 48 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00729-1, IRX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83800-4-PBS. Detection antibody: 83800-3-PBS. Standard: Ag19586. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00729-3, IRX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83800-2-PBS. Detection antibody: 83800-3-PBS. Standard: Ag19586. Range: 0.313-40 ng/mL