

For Research Use Only

IRX6 Recombinant antibody

Catalog Number:83800-5-RR



Basic Information

Catalog Number: 83800-5-RR	GenBank Accession Number: BC136573	Purification Method: Protein A purification
Size: 1000 ug/ml	GeneID (NCBI): 79190	CloneNo.: 240876D11
Source: Rabbit	UNIPROT ID: P78412	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG	Full Name: iroquois homeobox 6	
Immunogen Catalog Number: AG19586	Calculated MW: 446 aa, 48 kDa	
	Observed MW: 50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, NCCIT cells, HepG2 cells
Species Specificity: human	

Background Information

IRX6 is a member of the Iroquois homeobox gene family. Members of this family appear to play multiple roles during pattern formation of vertebrate embryos.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83800-5-RR (IRX6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83800-5-RR against Human IRX6 were performed. The affinity constant is 1.56 nM.